



**Chandler Loop example:
Reduced Thrombogenicity with Qvanteq's QSTH Surface
Technology applied onto Neurovascular Stents**

Kunze E., Rosenberger J.

December 20th, 2019
Arik Zucker, Stefano Buzzi

Experimental setup & method of dynamic human blood loop test

Test setup

- Flow loop model: Chandler loop
- Devices / run: Regular neurovascular stent vs. surface-treated stent per run and per donor
- Material: NiTi braiding; used for neurovascular stents

Test parameters


- Blood: human, freshly drawn (blood donation center Zurich; close proximity to Chandler Loop)
- Heparinization: 0.5 IU/ml
- Tube: Ø 4.75mm x 500mm
- Flow: 30 rpm
- Duration: 60 min
- Temp: 37°C

Analysis after runtime


- Blood extraction from loop tubes
 - Capturing **free-floating thrombi** (if present) with a filter
 - Determining **Thrombin-Antithrombin biomarker** to assess blood activation level difference between control & surface-treated stent, including also an empty tube as comparison
- Stent extraction from loop tubes rinsing in PBS and fixing adherent cells with PFA

Working principle of dynamic human blood loop test

Loop tube



Chandler loop setup



QSTH surface technology significantly reduces blood activation (TAT) & thrombi

DYNAMIC Human blood test in Chandler Loop	NiTi Stent (control)			NiTi Stent + QSTH surface			TAT ratio: control vs. surface treated	Blood loop without stent TAT [ug/ml]
	Stent after flow loop run	Free-floating thrombi in blood loop	TAT [ug/ml]	Stent after flow loop run	Free-floating thrombi in blood loop	TAT [ug/ml]		
	Donor CL77			22.1				
Donor CL79			16.2			0.06	270	0.06
Donor CL80			5.0			0.03	167	0.03
Avg ± Std			14.4	0.05			288	0.04

- QSTH technology
 - makes NiTi essentially “invisible” towards blood
 - mainly stent design influences thrombogenicity → hemodynamics



For further enquiries, please contact:

Arik Zucker, CEO

arik.zucker@qvanteq.com

Phone: +41 (0)44 500 98 90

Mobile: +41 (0)79 786 69 19

Qvanteq AG

Technoparkstrasse 1

CH – 8005 Zürich

Switzerland